

THE NEW PATIENT SPECIFIC BONE GRAFT

MyBone is a 3D printed porous implant made of a 100% biocompatible calcium-phosphate : Hydroxyapatite, that is the main component of the natural bone.


ADVANTAGES OF MYBONE



- PATIENT SPECIFIC IMPLANT
- HYDROXYAPATITE 100% BIOCOMPATIBLE
- STIMULATES BONE REGENERATION
- LOW INVASIVE SURGERY
- PATIENT'S COMFORT IMPROVED



3Dimensionnall tissue ingrow & regeneration
7 TIMES MORE EFFICIENT than the gold standard



THE 4 WEEK PROCESS OF MYBONE

- PATIENT'S CT-SCAN
- SURGEON'S REQUIREMENTS
- IMPLANT VALIDATION
- 3D PRINTING
- STERILIZATION
- READY TO USE

